


# VITA Pan Classic zur Bestimmung der natürlichen Zahnfarbe



# Mit Sonicare Healthy White um 2 Grade hellere Zähne in 2 Wochen!

**Conclusion** Sonicare HealthyWhite was shown to be effective in removing commonly observed extrinsic stain-forming pigments from tooth surfaces. An improvement of two Vitapan Classical shades was seen following 2, 3 and 6 weeks product use.



Before



After two weeks of use  
with HealthyWhite

Putt M, Milleman J, Jenkins W, Wei J, Schmitt P, Strate J. Data on file, 2007.